

CEL MXD

Independence Solution

The Tunstall CEL MXD is a Medical Alert System which utilizes a cellular network to communicate with our Emergency Response Center. A landline or traditional phone service is not required, nor does it require a cellular plan. The CEL MXD system includes a two-way voice base console and one personal help button which can be worn as a lanyard or wristband.

How do they work?

Subscribers call for help by pressing the button on their personal activator or CEL MXD console. The call is sent to a professional 24-hour Emergency Response Center where trained operators can access call details (subscriber's name, address, and relevant information) that are displayed on their screen.

The Emergency Response Center operator communicates with the subscriber through the CEL MXD console's loudspeaker and microphone to determine the nature of the call and the appropriate course of action.



Who are they for?

A wide variety of subscribers can benefit from the CEL MXD's state-of-the-art technology, including people who:

- Have been discharged from the hospital and require additional support and assistance at home
- Are of any age living alone
- Are older or have limited mobility
- Are subject to domestic violence or harassment
- Have household security issues or are living in fear of crime
- Are caring for someone in their home
- Have long-term conditions such as dementia, heart disease, or COPD.
- Have mental, physical, or sensory impairments

The CEL MXD offers family, caregivers, and friends the reassurance that their loved ones can use their personal activator to easily call for help whenever they need it in and around their home.

CEL MXD

Key Features

- Professional Emergency Response Centers
- Audible announcements
- 60 hour back-up console battery
- Integrated loud speaker & microphone
- Operates on the Rogers, Bell, Telus, or SAKSTEL cellular networks
- Supports up to 16 personal help buttons
- ETL listed
- Compatible with the MyActive Alert Fall Detector

Personal Help Button

- Multiple wearing options
- 600ft activator range
- Water resistant (IP67)

Technical Details

- Part(s) for this solution
- Dimensions
- Weight
- Battery
- Installation Requirements
- Maintenance Requirements

Independence Solution

Benefits

- Available 24/7, 365
- Subscriber is informed of the call progress and status
- Provides reassurance in the event of a power outage
- Enables two-way voice communication between subscriber and Emergency Call Response Operator
- No need for land-line, no cellular plan needed
- Allows multiple subscribers to use a single console unit
- Tunstall CEL MXD is FDA registered, complies with all medical device regulations, and is ETL listed to conform with UL 1635 and UL 1637 product safety standards.
- Automatically senses when a fall occurs

- Necklace or bracelet
- Range suitable for typical dwellings/homes
- Can be worn in shower or bath

- CEL MXD (MT100) Standard Activator:
(MT004) MyActive Fall Detector: (MT002)
- 2.5" x 5.25" x 7.8" (HxWxL)
- 1.8 lbs
- 60 hour back-up base battery
- AC power outlet and cellular coverage
- Periodic testing recommended