

myActive Alert

Independence Solution

Tunstall

SOLUTION SHEET

What is it?

The Tunstall myActive Alert features patented and clinically validated automatic fall detection algorithms. The myActive Alert also incorporates a water resistant HELP button to request help at the press of a button.

How does it work?

By pressing the button on the pendant, an emergency call is generated through the console unit and sent to the 24-hour Emergency Response Center (ERC). Trained operators can access call details (subscriber's name, address, and relevant information) that are displayed on their screen.

The ERC operator communicates with the subscriber through the console's loudspeaker and microphone to determine the nature of the call and the appropriate course of action.

Built-in automatic fall detection senses when the subscriber has fallen and triggers a call to the ERC.



Who is it for?

A wide variety of subscribers can benefit from the myActive personal emergency response pendant, including people who:

- Are at high risk for falls
- Have muscle weakness, especially in their legs
- Take medications that may cause dizziness
- Have hypotension (low blood pressure)
- Have neurological conditions such as Parkinson's Disease or Epilepsy
- Have mental, physical, or sensory impairments

The myActive offers family, caregivers, and friends the reassurance that their loved ones can use their myActive pendant to easily call for help whenever they need it from anywhere in and around the home.

The myActive Alert is compatible with the CEL MXD.

myActive Independence Solution

Key Features

- Advanced algorithms
- Auto cancel feature
- Small, lightweight design
- Compatibility
- 600ft radio range in free space
- Water resistant (IP67)
- Supervised auto presence

Benefits

- Enables the myActive to distinguish between fall and normal activities, such as sitting down in a chair
- Automatically cancels the call if steps are detected immediately after a fall event, eliminating unnecessary calls
- Discreet and easy to wear
- The myActive can be easily programmed to the CEL MXD
- Range suitable for in and around the home
- Can be worn in the shower or bath
- Notifies the Emergency Response Center if the myActive is not detected within range of the console for a specified period of time

Specifications

• Part Number:	MT002
• Dimension:	2.165" x 1.468" x 0.603"
• Weight:	1oz
• RF Transmission Frequency:	418 MHz
• Operating Temperature:	-10°C to 60°C
• Battery Type:	CR2450 Lithium
• Battery Life:	12-18 months typical
• FCC ID:	NY5MAA01
• IC:	9663A-MAA01
• Product Compatibility:	CEL MXD

For additional information please call
800.892.2205
or email info.canada@tunstall.com
111 Zenway Blvd, Unit #6A
Woodbridge, Ontario L4H3H9 Canada
canada.tunstall.com

Tunstall

Due to the many types of falls that can occur, the user should always press the HELP button whenever help is needed and not rely solely upon the determination of the fall detection