

# iVi

## Intelligent Pendant with Fall Detector

The Fall Detector is a unique personal device. It provides reliable detection of a potentially serious fall by the wearer and then triggers a call for help via Caresse+ or Vi+ personal emergency response unit. The device is attractively styled while being small and light enough to be worn unobtrusively.



### Features

- Attractive and discreet
- Ergonomically styled with matte finish
- Manual trigger button
- Activation indicator light and audible alarm or sound
- User replaceable battery
- Belt clip or holster wearing options
- Easy to keep clean
- 4 hour automatic transmission of status reports for pendant supervision
- Wireless connection to Caresse+ and Vi+
- Generates an audible warning as well as triggering a call on the host unit

# Technical Details

Part Number	P66004/07 (light grey)
Battery	Supplied in a box with user manual, battery and additional spare battery, belt and holster, and table stand.
Radio Frequency	312.000MHz
Battery	<ul style="list-style-type: none"> <li>• CR2450 (3v Lithium cell). Only use CR2450 batteries manufactured by Sony. Typical life up to 12 months.</li> <li>• When the battery is nearing the end of its life, the light flashes when the button is pressed and the call recipient is automatically notified of the low battery state. Battery status is also reported automatically during the 4 hourly supervision transmissions.</li> </ul>
Dimensions	58x38x14mm [2.28"x1.5"x0.55"]
Weight	25g [0.88oz]
Operating Temp Range	(to perform to full specification) = 10° C -45° C, storage = -10° C to 50° C
Radio Range	Up to 120m [390ft] (line of sight)
Supervision	Auto Low Battery
Standards	UL1637, FCC CFR 47, Part 15 FDA Class II Medical Device , Canadian RSS-102
Product Compatibility	Caresse+ and Vi+

Tunstall Canada Inc.  
111 Zenway Blvd, Unit 6A  
Woodbridge, Ontario  
Canada L4H 3H9

Tel: 1-800-892-2205  
Tel: 905-677-1144  
Email: [info.canada@tunstall.com](mailto:info.canada@tunstall.com)  
[canada.tunstall.com](http://canada.tunstall.com)  
[americas.tunstall.com](http://americas.tunstall.com)